

## Reaction vessels, spherical, with thermostatic jacket and withdrawal valve



Reaction vessels, spherical, with thermostatic jacket and withdrawal valve

Without graduation, with laboratory flange (DN), with groove, with two connectors (with threads GL 14 / GL 18 or ball flanges KF 15), made of DURAN® tubing, other sizes on request, delivery with 2 plastic olives GL14/18 and screw caps

Nom. size (DN)	Capacity (l)	Vessel dia. Outer (mm)	Vessel dia. Inner (mm)	Total height (approx. mm)	Connectors	Valve bore (mm)	Pack (qty.)	Cat. No.
60	0,25	105	85	255	GL 14	10	1	6 0401 49
60	0,5	131	105	280	GL 14	10	1	6 0401 58
60	1	165	131	315	GL 14	10	1	6 0401 70
60	0,25	105	85	255	KF 15	10	1	6 0411 49
60	0,5	131	105	280	KF 15	10	1	6 0411 58
60	1	165	131	315	KF 15	10	1	6 0411 70
100	2	206	165	360	GL 18	10	1	6 0402 73
100	4	236	206	395	GL 18	10	1	6 0402 79

100	6	280	236	440	GL 18	10	1	6 0402 85
100	2	206	165	360	KF 15	10	1	6 0412 73
100	4	236	206	395	KF 15	10	1	6 0412 79
100	6	280	236	440	KF 15	10	1	6 0412 85
150	6	280	236	455	GL 18	10	1	6 0404 85
150	6	280	236	455	KF 15	10	1	6 0414 85
150	10	350	275	525	KF 15	20	1	6 0414 88
200	6	280	236	440	KF 15	10	1	6 0415 85
200	10	350	275	525	KF 15	20	1	6 0415 88

Generated on 2026-05-18