

## Round bottom flasks, wide neck



Round bottom flasks, wide neck

Acc. to DIN ISO 24450 and in addition to standard, with beaded rim, material: DURAN®

Capacity (ml)	Flask dia. (mm)	Neck dia. (mm)	Height (mm)	Acc. to DIN	Pack (qty.)	Cat. No.
50	51	34	105	-	10	7 0350 28
100	64	35	110	DIN ISO 24450	10	7 0350 37
250	85	51	143	DIN ISO 24450	10	7 0350 49
500	105	50	168	DIN ISO 24450	10	7 0350 58
1000	131	50	200	DIN ISO 24450	10	7 0350 70
1000	131	65	200	-	10	7 0350 71
2000	165	76	240	DIN ISO 24450	10	7 0350 72
2000	166	50	240	-	10	7 0350 73
3000	185	65	260	-	1	7 0350 76
4000	206	76	290	DIN ISO 24450	1	7 0350 79
5000	223	65	310	-	1	7 0350 82
6000	236	89	330	DIN ISO 24450	1	7 0350 84
6000	236	65	330	-	1	7 0350 85
10000	279	89	420	-	1	7 0350 88
20000	345	89	520	-	1	7 0350 94