

Round-bottom flasks with conical ground joint, amber



Round-bottom flasks with conical ground joint, amber

According to DIN EN ISO 4797, with inscription label, material: DURAN®

Capacity (ml)	Socket (NS)	Flask outer dia. (mm)	Pack (qty.)	Cat. No.
50	29/32	51	10	3 0031 28
100	29/32	64	10	3 0031 37
250	29/32	85	10	3 0031 49
500	29/32	105	10	3 0031 58
1000	29/32	131	10	3 0031 70