

## Screwthread tubes with NS cone



Screwthread tubes with NS cone

Acc. to DIN 12257 and in addition to standard, with DIN thread (GL) and cone (NS), made of DURAN® tubing

Cone (NS)	Thread (GL)	Acc. to DIN	Pack (quantity)	Cat. No.
14/23	14	DIN 12257	1	1 3114 14
19/26	14	DIN 12257	1	1 3114 19
24/29	14	-	1	1 3114 24
29/32	14	-	1	1 3114 29
14/23	18	-	1	1 3118 14
19/26	18	DIN 12257	1	1 3118 19
24/29	18	-	1	1 3118 24
29/32	18	DIN 12257	1	1 3118 29
14/23	25	-	1	1 3125 14
19/26	25	-	1	1 3125 19
24/29	25	-	1	1 3125 24
29/32	25	DIN 12257	1	1 3125 29
29/32	32	DIN 12257	1	1 3132 29