

## Sedimentation cone acc. to Imhoff, without stopcock



Sedimentation cone acc. to Imhoff, without stopcock

Sedimentation cone acc. to DIN 12672,  
Scale: up to 100 ml,  
ring mark: at 1000 ml,  
subdivision: 0,1 ml (0-2 ml), 0,5 ml (2-10 ml),  
1 ml (10-40 ml), 2 ml (20-100 ml),  
material: borosilicate glass 3.3  
Cat. No.: 5 6700 01 (pack: 1 piece)

