

Sedimentation cone acc. to Imhoff, without stopcock



Sedimentation cone acc. to Imhoff, without stopcock

Sedimentation cone acc. to DIN 12672, Scale: up to 1000 ml, ring mark: at 1000 ml, subdivision: 0,1 ml (0-2 ml), 0,5 ml (2-10 ml), 1 ml (10-40 ml), 2 ml (20-100 ml), 50 ml (100-1000 ml), material: borosilicate glass 3.3
Cat. No.: 5 6700 03 (pack: 1 piece)

