

Separating funnels, acc. to Squibb, amber, with glass plug, without scale



Separating funnels, acc. to Squibb, amber, with glass plug, without scale

With NS socket and PE stopper, without scale, with solid clear glass plugs, with screwthread retaining device, material: borosilicate glass 3.3

Capacity (ml)	Socket (NS)	Bore (mm)	Pack (qty.)	Cat. No.
50	19/26	2,5	2	4 0180 28
100	19/26	2,5	2	4 0180 37
250	29/32	4,0	2	4 0180 49
500	29/32	4,0	2	4 0180 58
1000	29/32	6,0	2	4 0180 70
2000	29/32	6,0	1	4 0180 73