

Separating funnels, acc. to Squibb, with glass plug, with scale



Separating funnels, acc. to Squibb, with glass plug, with scale

With NS socket and PE stopper, with solid glass plug, with screwthread retaining device, with scale (approx. volume / ml), material: borosilicate glass 3.3

| Capacity (ml) | Socket (NS) | Bore (mm) | Scale (ml) | Pack (qty.) | Cat. No. |
|---------------|-------------|-----------|------------|-------------|-----------|
| 50 | 19/26 | 2,5 | 2 | 2 | 4 0120 28 |
| 100 | 19/26 | 2,5 | 5 | 2 | 4 0120 37 |
| 250 | 29/32 | 4,0 | 5 | 2 | 4 0120 49 |
| 500 | 29/32 | 4,0 | 10 | 2 | 4 0120 58 |
| 1000 | 29/32 | 6,0 | 20 | 2 | 4 0120 70 |
| 2000 | 29/32 | 6,0 | 50 | 1 | 4 0120 73 |