

## Socket members



Socket members

Acc. to DIN 12244 and in addition to standard, made of DURAN® tubing

Socket (S)	Total length approx. (mm)	Acc. to DIN	Pack (quantity)	Cat. No.
13/2	100	-	10	1 2013 02
13/5	100	DIN 12244	10	1 2013 05
19	120	DIN 12244	10	1 2019 09
29	120	DIN 12244	10	1 2029 15
35	120	DIN 12244	10	1 2035 20
40	120	DIN 12244	1	1 2040 25
41	120	DIN 12244	1	1 2041 25
51	150	DIN 12244	1	1 2051 30
65	170	DIN 12244	1	1 2065 40