

## Sockets with special tube dimensions



Sockets with special tube dimensions

Conical ground sockets, made of DURAN® tubing

Socket (NS)	Dia. (mm)	Wall thickn. (mm)	Total length approx. (mm)	Pack (qty.)	Cat. No.
14/23	16	1,2	160	10	1 1010 01
14/23	22	1,8	145	10	1 1010 02
29/32	33	2,0	210	10	1 1010 05
29/32	40	2,3	380	10	1 1010 07
45/50	48	2,3	280	10	1 1010 11