



Spherical Joints

Our spherical joints are manufactured to inch dimensions.
The measures are in accordance
with DIN 12244 and 12264.

{loadposition user3}

Nominal size	d1 (mm)	d2 (mm)	d3 (mm)	d4 (mm)
S13	12,700	12,706	12,50	12,5
S19	19,050	19,057	18,90	18,9
S29	28,575	28,582	28,35	28,3
S35	34,925	34,934	34,65	34,6
S40	38,100	38,109	37,80	37,7
S41	41,275	41,284	40,95	40,9
S51	50,800	50,810	49,40	49,3
S64	63,500	63,510	62,95	62,9

Generated on 2026-05-25