

Sprayer



Sprayer

Complete sprayer consisting of: Erlenmeyer flask, (Capacity: 100 ml), with socket NS 19/26, material: DURAN®, sprayer head (made of DURAN® tubing), with clamp and rubber bulb, packaging unit: 1 piece Cat. No.: 5 4700 00

Components Description	Pack (qty.)	Cat. No.
Sprayer head	1	5 4700 04
Flask, 100 ml, NS 19/26	1	3 0219 37
Clamp NS 19/26	1	1 1037 19
Rubber bulb	1	5 4700 05