

Stirrer bearings with dynamic double seal



Stirrer bearings with dynamic double seal

Stirrer bearings for applications under vacuum or under slight overpressure. With extended guidance for reduced vibrations and a longer lifetime. Material: borosilicate glass 3.3 (closure), PTFE (sealings)

| For shaft (dia. mm) | Joint (NS) | Pack (qty.) | Cat. No. |
|---------------------|------------|-------------|-----------|
| 10 | 29/32 | 1 | 6 3060 04 |
| 10 | 45/40 | 1 | 6 3060 06 |
| 16 | 29/32 | 1 | 6 3062 04 |
| 16 | 45/40 | 1 | 6 3062 06 |

Other joints on request