

Stirrer with stirrer bearing



Stirrer with stirrer bearing

Shaft with blade and PTFE security screw, gastight closure of vacuum stirring apparatuses up to 10^{-3} mbar at a continuous operation speed of max. 800 rpm, for round bottom flasks 250-2000 ml, shaft made of DURAN® tubing

Description	Pack (qty.)	Cat. No.
Stirrer with stirrer bearing, complete	1	6 3100 00
Shaft, Dia. 10 x L 400 mm	1	6 3105 01
Blade, PTFE, 70 x 24 x 3 mm	1	6 3105 02
Security screw, PTFE (see picture)	1	6 3105 03
Closure, thread GL 25, cone NS 29/32	1	6 3105 04
Sealing, PTFE/graphite, 1 set	1	6 3105 05
Pressure spring, stainless steel	1	6 3105 06
Screw cap GL 25, PP	1	6 3105 07