

Stirrers with lateral blade and GL thread connection



Stirrers with lateral blade and GL thread connection

Made of precision calibrated glass tube, length of the ground and polished bearing surface: 160 mm, blade and security screw: PTFE, shaft made of DURAN® tubing

Spare parts

For joint (NS)	Shaft dia. x L (mm)	Blade (mm)	Thread (GL)	Pack (qty.)	Cat. No.
29/32	10 x 310	50 x 18 x 3	10	1	6 3070 31
29/32	10 x 310	68 x 24 x 3	10	1	6 3070 33
29/32	10 x 390	75 x 24 x 3	10	1	6 3070 35
45/40	16 x 460	75 x 24 x 3	14	1	6 3070 68
29/32	10 x 480	125 x 24 x 3	10	1	6 3070 30
45/40	16 x 520	125 x 35 x 3	14	1	6 3070 61
45/40	16 x 560	125 x 35 x 3	14	1	6 3070 63
45/40	16 x 620	125 x 35 x 3	14	1	6 3070 65

Description	Pack (qty.)	Cat. No.
Shaft 10 x 310 with journal 8 x 40	1	6 3072 01
Shaft 10 x 390 with journal 8 x 40	1	6 3072 05
Shaft 10 x 480 with journal 8 x 40	1	6 3072 10
Shaft 16 x 460 with journal 14 x 50	1	6 3072 08

Shaft 16 x 520 with journal 14 x 50	1	6 3072 11
Shaft 16 x 560 with journal 14 x 50	1	6 3072 13
Shaft 16 x 620 with journal 14 x 50	1	6 3072 15
Shaft 10 x 310 with thread GL 10	1	6 3072 51
Shaft 10 x 390 with thread GL 10	1	6 3072 55
Shaft 10 x 480 with thread GL 10	1	6 3072 68
Shaft 16 x 460 with thread GL 14	1	6 3072 60
Shaft 16 x 520 with thread GL 14	1	6 3072 71
Shaft 16 x 560 with thread GL 14	1	6 3072 73
Shaft 16 x 620 with thread GL 14	1	6 3072 75
Blade 50 x 18 x 3, PTFE	1	6 3073 01
Blade 68 x 24 x 3, PTFE	1	6 3073 03
Blade 75 x 24 x 3, PTFE	1	6 3073 05
Blade 125 x 24 x 3, PTFE	1	6 3073 10
Blade 125 x 35 x 3, PTFE	1	6 3073 11
PTFE security screw	1	6 3105 03

Generated on 2026-05-18