

Stirrers with lateral blade and journal



Stirrers with lateral blade and journal

Made of precision calibrated glass tube, length of the grounded and polished bearing surface: 160 mm, blade and security screw: PTFE, shaft made of DURAN® tubing

Spare parts

| For joint (NS) | Shaft dia. x L (mm) | Blade (mm) | Journal dia. x L (mm) | Pack (qty.) | Cat. No. |
|----------------|---------------------|--------------|-----------------------|-------------|-----------|
| 29/32 | 10 x 310 | 50 x 18 x 3 | 8 x 40 | 1 | 6 3070 01 |
| 29/32 | 10 x 310 | 68 x 24 x 3 | 8 x 40 | 1 | 6 3070 03 |
| 29/32 | 10 x 390 | 75 x 24 x 3 | 8 x 40 | 1 | 6 3070 05 |
| 45/40 | 16 x 460 | 75 x 24 x 3 | 10 x 50 | 1 | 6 3070 08 |
| 29/32 | 10 x 480 | 125 x 24 x 3 | 8 x 40 | 1 | 6 3070 10 |
| 45/40 | 16 x 520 | 125 x 35 x 3 | 10 x 50 | 1 | 6 3070 11 |
| 45/40 | 16 x 560 | 125 x 35 x 3 | 10 x 50 | 1 | 6 3070 13 |
| 45/40 | 16 x 620 | 125 x 35 x 3 | 10 x 50 | 1 | 6 3070 15 |

| Description | Pack (qty.) | Cat. No. |
|-------------------------------------|-------------|-----------|
| Shaft 10 x 310 with journal 8 x 40 | 1 | 6 3072 01 |
| Shaft 10 x 390 with journal 8 x 40 | 1 | 6 3072 05 |
| Shaft 10 x 480 with journal 8 x 40 | 1 | 6 3072 10 |
| Shaft 16 x 460 with journal 10 x 50 | 1 | 6 3072 08 |
| Shaft 16 x 520 with journal 10 x 50 | 1 | 6 3072 11 |
| Shaft 16 x 560 with journal 10 x 50 | 1 | 6 3072 13 |

| | | |
|-------------------------------------|---|-----------|
| Shaft 16 x 620 with journal 10 x 50 | 1 | 6 3072 15 |
| Shaft 10 x 310 with thread GL 10 | 1 | 6 3072 51 |
| Shaft 10 x 390 with thread GL 10 | 1 | 6 3072 55 |
| Shaft 10 x 480 with thread GL 10 | 1 | 6 3072 68 |
| Shaft 16 x 460 with thread GL 14 | 1 | 6 3072 60 |
| Shaft 16 x 520 with thread GL 14 | 1 | 6 3072 71 |
| Shaft 16 x 560 with thread GL 14 | 1 | 6 3072 73 |
| Shaft 16 x 620 with thread GL 14 | 1 | 6 3072 75 |
| Blade 50 x 18 x 3, PTFE | 1 | 6 3073 01 |
| Blade 68 x 24 x 3, PTFE | 1 | 6 3073 03 |
| Blade 75 x 24 x 3, PTFE | 1 | 6 3073 05 |
| Blade 125 x 24 x 3, PTFE | 1 | 6 3073 10 |
| Blade 125 x 35 x 3, PTFE | 1 | 6 3073 11 |
| PTFE security screw | 1 | 6 3074 00 |

Generated on 2026-05-18