


Stirrers with paddles and GL thread connection



Stirrers with paddles and GL thread connection

Made of precision calibrated glass tube, ground and polished bearing surface, shaft dia.: 10 mm, thread connection: GL 10, paddles and pin with disc: PTFE, shaft: made of DURAN® tubing

PTFE paddles

For joint (NS)	Shaft dia. Å? x L (mm)	Circle (mm)	Pack (Stck.)	Cat. No.
29/32	10 x 310	50	1	6 3080 31
29/32	10 x 310	90	1	6 3080 33
29/32	10 x 400	90	1	6 3080 35
29/32	10 x 460	90	1	6 3080 37
29/32	10 x 520	90	1	6 3080 39
45/40	10 x 560	120	1	6 3080 40
45/40	10 x 620	120	1	6 3080 42

Circle (mm)	Pack (qty.)	Cat. No.
50	1	6 3081 01
90	1	6 3081 03
120	1	6 3081 05