

Test tube sprayers



Test tube sprayers

Complete sprayer consisting of: Test tube and sprayer head (made of DURAN® tubing), with clamp and rubber bulb, packaging unit: 1 piece 6 ml, NS 14/23 Cat. No.: 5 4706 00 12 ml, NS 19/26 Cat. No.: 5 4712 00

Components Description	Pack (qty.)	Cat. No.
Sprayer head, 6 ml, NS 14/23	1	5 4706 05
Sprayer head, 12 ml, NS 19/26	1	5 4712 05
Test tube for head 6 ml, NS 14/23	1	5 4706 06
Test tube for head 12 ml, NS 19/26	1	5 4712 06
Clamp NS 14/23	1	1 1037 14
Clamp NS 19/26	1	1 1037 19
Rubber bulb	1	5 4700 05