

## Test tubes with NS joint, with graduation



Test tubes with NS joint, with graduation

Round bottom glass with NS socket, with PE stopper, with scale in high-contrast enamel, made of DURAN® tubing

Socket (NS)	Dia. x height (mm)	Capacity (ml)	Subdivision (ml)	Pack (qty.)	Cat. No.
12/21	15 x 165	10	0,1	10	5 7060 01
14/23	17 x 200	20	0,2	10	5 7060 02
14/23	17 x 220	25	0,2	10	5 7060 03
14/23	19 x 220	30	0,5	10	5 7060 04