

## Test tubes with NS joint without graduation



Test tubes with NS joint without graduation

Round bottom glass with NS socket, available without stopper, with PE stopper or with hexagonal hollow glass stopper, made of DURAN® tubing

Socket (NS)	Dia. x height (mm)	Description	Pack (qty.)	Cat. No.
14/23	16 x 100	Without stopper	10	5 7000 07
14/23	16 x 125	Without stopper	10	5 7000 10
14/23	16 x 160	Without stopper	10	5 7000 14
14/23	18 x 150	Without stopper	10	5 7000 26
14/23	18 x 180	Without stopper	10	5 7000 30
19/26	22 x 150	Without stopper	10	5 7000 34
24/29	28 x 150	Without stopper	10	5 7000 38
24/29	28 x 200	Without stopper	10	5 7000 42
14/23	16 x 100	PE stopper	10	5 7001 07
14/23	16 x 125	PE stopper	10	5 7001 10
14/23	16 x 160	PE stopper	10	5 7001 14
14/23	18 x 150	PE stopper	10	5 7001 26
14/23	18 x 180	PE stopper	10	5 7001 30
19/26	22 x 150	PE stopper	10	5 7001 34
24/29	28 x 150	PE stopper	10	5 7001 38
24/29	28 x 200	PE stopper	10	5 7001 42

14/23	16 x 100	Hexagonal glass stopper	10	5 7002 07
14/23	16 x 125	Hexagonal glass stopper	10	5 7002 10
14/23	16 x 160	Hexagonal glass stopper	10	5 7002 14
14/23	18 x 150	Hexagonal glass stopper	10	5 7002 26
14/23	18 x 180	Hexagonal glass stopper	10	5 7002 30
19/26	22 x 150	Hexagonal glass stopper	10	5 7002 34
24/29	28 x 150	Hexagonal glass stopper	10	5 7002 38
24/29	28 x 200	Hexagonal glass stopper	10	5 7002 42

Generated on 2026-05-19