

Three-neck round-bottom flasks, side necks parallel



Three-neck round-bottom flasks, side necks parallel

Reaction flasks acc. to DIN 12392 and in addition to DIN, with two parallel >side necks, with conical ground joints, with inscription label, material: DURAN®

Capacity (ml)	Center neck (NS)	Side neck (NS)	Flask outer dia. (mm)	Acc. to DIN	Pack (quantity)	Cat. No.
100	24/29	19/26	64	-	1	3 3032 37
250	24/29	19/26	85	-	1	3 3032 49
500	24/29	19/26	105	-	1	3 3032 58
1000	24/29	19/26	131	-	1	3 3032 70
2000	24/29	19/26	166	-	1	3 3032 73
100	29/32	14/23	64	-	1	3 3041 37
250	29/32	14/23	85	-	1	3 3041 49
500	29/32	14/23	105	-	1	3 3041 58
1000	29/32	14/23	131	-	1	3 3041 70
2000	29/32	14/23	166	-	1	3 3041 73
3000	29/32	14/23	185	-	1	3 3041 76
4000	29/32	14/23	207	-	1	3 3041 79
6000	29/32	14/23	236	-	1	3 3041 85
100	29/32	19/26	64	-	1	3 3042 37
250	29/32	19/26	85	-	1	3 3042 49

500	29/32	19/26	105	-	1	3 3042 58
1000	29/32	19/26	131	-	1	3 3042 70
2000	29/32	19/26	166	-	1	3 3042 73
100	29/32	29/32	64	-	1	3 3044 37
250	29/32	29/32	85	-	1	3 3044 49
500	29/32	29/32	105	-	1	3 3044 58
1000	29/32	29/32	131	DIN 12392	1	3 3044 70
2000	29/32	29/32	166	DIN 12392	1	3 3044 73
3000	29/32	29/32	185	-	1	3 3044 76
4000	29/32	29/32	207	-	1	3 3044 79
5000	29/32	29/32	223	-	1	3 3044 82
6000	29/32	29/32	236	-	1	3 3044 85
10000	29/32	29/32	279	-	1	3 3044 88
1000	45/40	29/32	131	-	1	3 3054 70
2000	45/40	29/32	166	DIN 12392	1	3 3054 73
3000	45/40	29/32	185	-	1	3 3054 76
4000	45/40	29/32	207	DIN 12392	1	3 3054 79
5000	45/40	29/32	223	-	1	3 3054 82
6000	45/40	29/32	236	DIN 12392	1	3 3054 85
10000	45/40	29/32	279	-	1	3 3054 88
20000	45/40	29/32	345	-	1	3 3054 94

Generated on 2026-05-18