

Tubes with socket and cone



Tubes with socket and cone

Conical ground sockets and cones, made of DURAN® tubing

| Socket (NS) | Cone (NS) | Dia. (mm) | Wall thicken. (mm) | Total length approx. (mm) | Pack (qty.) | Cat. No. |
|-------------|-----------|-----------|--------------------|---------------------------|-------------|-----------|
| 14/23 | 14/23 | 17 | 1,8 | 120 | 1 | 1 1014 11 |
| 14/23 | 29/32 | 26 | 2,0 | 130 | 1 | 1 1014 14 |
| 19/26 | 19/26 | 22 | 1,8 | 125 | 1 | 1 1014 22 |
| 29/32 | 29/32 | 33 | 2,0 | 150 | 1 | 1 1014 44 |