

Two-neck pear-shaped flasks



Two-neck pear-shaped flasks

Acc. to DIN 12383 and in addition to DIN standard, with angled side neck (25°), with conical ground joints, with inscription label, material: DURAN®

Capacity (ml)	Center neck (NS)	Side neck (NS)	Flask outer dia. (mm)	Acc. to DIN	Pack (quantity)	Cat. No.
25	14/23	14/23	40	-	1	3 2211 22
50	14/23	14/23	50	-	1	3 2211 28
100	14/23	14/23	64	-	1	3 2211 37
250	14/23	14/23	85	-	1	3 2211 49
25	19/26	14/23	40	-	1	3 2221 22
50	19/26	14/23	50	DIN 12383	1	3 2221 28
100	19/26	14/23	64	DIN 12383	1	3 2221 37
250	19/26	14/23	85	DIN 12383	1	3 2221 49