

Vacuum stopcocks, 3-way



Vacuum stopcocks, 3-way

Acc. to DIN 12554 and in addition to DIN standard, T-shape, with right-angled side arms and hollow glass plug, material: borosilicate glass 3.3, leakage rate acc. to DIN 12540-1

Stopcock(NS)	Bore (mm)	Acc. to DIN	Pack (qty.)	Cat. No.
14,5	2,5	DIN 12545	1	2 9614 02
18,8	4,0	DIN 12545	1	2 9618 04
21,5	6,0	DIN 12545	1	2 9621 06
24	8,0	(1)	1	2 9624 08
29	10,0	DIN 12554	1	2 9629 10

(1) in addition to DIN standard