

Water estimators acc. to Dean Stark



Water estimators acc. to Dean Stark

Water estimators acc. to DIN 12420 (xylool method), capacity: 10:0,1 ml, with round bottom flask 500 ml, ground joint connections: NS29/32, measuring tube with or without stopcock, with Liebig condenser, jacket length: 400 mm, material: measuring tube and condenser made of DURAN® tubing, flask: DURAN®

Zero point (position)	Stopcock	Pack (qty.)	Cat. No.
at the stopcock	with	1	5 6000 00
above stopcock	with	1	5 6000 02
at the tip	without	1	5 6000 04

Components	Pack (qty.)	Cat. No.
Measuring tube with stopcock (zero point at stopcock)	1	5 6003 01
Meas. tube with stopcock (zero point above stopcock)	1	5 6003 03
Meas. tube without stopcock (zero point at the tip)	1	5 6003 05
Round bottom flask, 500 ml, NS 29/32	1	3 0029 58
Liebig condenser, jacket length: 400 mm, NS 29/32	1	5 2104 04