

Weighing bottles with NS stopper



Weighing bottles with NS stopper

With interchangeable NS stopper, with inscription label, with fused-in size (mm) and capacity (ml), available in tall or low form, stopper: DURAN®, bottle made of DURAN® tubing

Dimensions Dia. x H (mm)	Joint (NS)	Shape	Capacity (ml)	Pack (qty.)	Cat. No.
20 x 40	19/12	tall	5	6	5 5600 02
25 x 40	24/12	tall	10	6	5 5600 05
30 x 50	29/12	tall	20	6	5 5600 08
35 x 50	34/12	tall	30	6	5 5600 12
40 x 50	40/12	tall	40	6	5 5600 15
40 x 65	40/12	tall	45	6	5 5600 19
35 x 70	34/12	tall	45	6	5 5600 22
40 x 80	40/12	tall	70	6	5 5600 25
50 x 80	50/12	tall	110	10	5 5600 28
25 x 25	24/12	low	6	6	5 5605 02
40 x 25	40/12	low	15	6	5 5605 04
35 x 30	34/12	low	15	6	5 5605 09
50 x 30	50/12	low	30	6	5 5605 13
60 x 30	60/12	low	45	6	5 5605 16
80 x 30	80/12	low	80	6	5 5605 21