

## Übergangs-/Expansionsstücke Hülse/Kern



### Übergangs-/Expansionsstücke Hülse/Kern

Nach DIN 12257 und in Ergänzung zur DIN-Norm, geschliffene Kegelhülsen und Kegelkerne, hergestellt aus DURAN®-Rohr

Hülse (NS)	Kern (NS)	Nach DIN	VPE (Stück)	Art.-Nr.
14/23	12/21	-	1	1 5014 12
19/26	14/23	DIN 12257	1	1 5019 14
24/29	14/23	-	1	1 5024 14
24/29	19/26	-	1	1 5024 19
29/32	14/23	DIN 12257	1	1 5029 14
29/32	19/26	DIN 12257	1	1 5029 19
29/32	24/29	-	1	1 5029 24
34/35	19/26	-	1	1 5034 19
34/35	24/29	-	1	1 5034 24
34/35	29/32	-	1	1 5034 29
45/40	29/32	-	1	1 5045 29
45/40	34/35	-	1	1 5045 34
60/46	29/32	-	1	1 5060 29
60/46	45/40	-	1	1 5060 45
71/51	29/32	-	1	1 5071 29
71/51	45/40	-	1	1 5071 45

