

Übergangs-/Reduzierstücke Hülse/Kern



Übergangs-/Reduzierstücke Hülse/Kern

Nach DIN 12257 und in Ergänzung zur DIN-Norm, geschliffene Kegelhülsen und Kegelkerne, hergestellt aus DURAN®-Rohr

Hülse (NS)	Kern (NS)	Nach DIN	VPE (Stück)	Art.-Nr.
12/21	14/23	-	1	1 5012 14
14/23	19/26	DIN 12257	1	1 5014 19
14/23	24/29	-	1	1 5014 24
14/23	29/32	DIN 12257	1	1 5014 29
19/26	24/29	-	1	1 5019 24
19/26	29/32	DIN 12257	1	1 5019 29
19/26	34/35	-	1	1 5019 34
24/29	29/32	-	1	1 5024 29
24/29	34/35	-	1	1 5024 34
24/29	45/40	-	1	1 5024 45
29/32	34/35	-	1	1 5029 34
29/32	45/40	DIN 12257	1	1 5029 45
29/32	60/46	-	1	1 5029 60
29/32	71/51	-	1	1 5029 71
45/40	60/46	-	1	1 5045 60
45/40	71/51	-	1	1 5045 71

